



PATIENT
SUPPORT
SYSTEMS

MEDICAL • DESIGN • INNOVATION



Theraflow™
Alternating Mattress System

Comfort focused pressure therapy



Supports tissue reperfusion
with 1-in-3 alternating cycle.



Adapts support automatically
with intelligent weight detection.



Protects vulnerable pressure areas
with elevated heel zone.




Supports infection control
with fully welded construction.



Key features :

- ▶ Supportive 1-in-3 alternating pressure therapy
- ▶ Hyper-Stretch comfort coverlet included
- ▶ Quiet 10 LPM intelligent control unit
- ▶ Elevated heel pressure relief
- ▶ Fully welded infection control construction



Top cover & base removed for illustration purposes



PATIENT
SUPPORT
SYSTEMS

MEDICAL • DESIGN • INNOVATION

Unit C2, 15 Narabang Way Belrose, Sydney NSW 2085

☎ 1800 432 247 🖨 02 9844 5445 ✉ service@patientsupportsystems.com

www.patientsupportsystems.com

TheraflowTM

Alternating Mattress System

Specifications

Control Unit:

Dimensions:	W 27 × H 21 × D 10 cm
Weight	2.5kg
Mode of Operation	Non-Continuous
Power Rating	240V 50Hz 12VA
Transport Function	Hose End Cap
Air Flow	10 Litres per Minute
Auto Startup	✓
Auto Weight Detect	✓
Digital Comfort Control	✓
Self Diagnostics	✓
Audible/Visual alarm	✓
Mute Function	✓
Static Mode	Timeout - 60 minutes
Care Mode	Timeout - 20 minutes
Manual Incline Mode	✓
Auto Anti Tamper-control panel lock function	✓
CE Certification and C-Tick Listed N27882	✓

Warranty

- ▶ Control Unit: 2 years
- ▶ Soft Goods: 2 years

ARTG: 175810



Mattress Replacement :

Dimensions Standard	L 200 x W 90 x H 20 cm
Dimensions King	L 200 x W 107 x H 20 cm
Alternating Cell Type	Single Cell
Overall Cells Count	20
Cell Material	Pure TPU
Base/s	Assorted PU and PVCs, welded
Cell Cycle	1 in 3
Alternating Mode	10 minutes per alternation
Infection Control	Welded Top Cover, Zip and Base
Mattress Attachment	Non Slip PVC Straps

Hyper-stretch coverlet :

Dimensions	L 200 x W 90 x H 1cm
Attachment	Zips to PVC flaps around Cell Section
Cover Material	High Performance Micromesh
Order code	TF-08-1H-LS-090-NE7-CMP